

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
 ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321
 13088 Berlin, Germany

T: +49 30 927 030-0
 F: +49 30 927 030-30
 info@sifin.de
 www.sifin.de

| | |
|--------------------|--|
| Product | Anti-Salmonella H:6 |
| Catalogue no. | TR 1435 |
| Package size | 1 ml |
| Lot no. | 2820318 |
| Expiry date | 2021-03 |
| Intended use | Anti-Salmonella H:6 records all <i>Salmonella</i> strains, as long as they contain antigen H:1,6 in slide agglutination. |
| Specificity | No agglutination with heterologous <i>Salmonella</i> H antigens and <i>Salmonella</i> O antigens. |
| Preservatives | Sodium azide 0.9 mg/ml |
| Storage | 2 °C...8 °C |

Berlin,

 place, date

23.05.2018



 signature Head of Quality Control