

Certificate of Quality Control

Manufacturer **SIFIN** Institut für Immunpräparate und Nährmedien GmbH Berlin,
Berliner Allee 317-321, 13088 Berlin, Germany

Product **Enteroclon Anti-Salmonella H 1**

Catalogue No.
(Package size) TR 1437 (1 ml)

Lot No. **180 11 12**

Expiry date **2015-11**

Intended use Enteroclon Anti-Salmonella H 1 records all Salmonella strains, as long as they contain antigen H 1,2, H 1,5, H 1,6, H 1,7.

Specificity No agglutination with heterologous Salmonella H antigens and Salmonella O antigens.

Preservatives Sodium azide 0.9 mg/ml

Storage 2...8 °C

Berlin, 04.01.13

K. P. G. W.
Head of Quality Control