

## Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN  
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485


sifin diagnostics gmbh

Berliner Allee 317-321  
13088 Berlin, Germany

T: +49 30 927 030-0  
F: +49 30 927 030-30  
info@sifin.de  
www.sifin.de

<b>Product</b>	<b>Anti-Salmonella H:1</b>
Catalogue no.	TR 1437
Package size	1 ml
<b>Lot no.</b>	<b>187 10 14</b>
<b>Expiry date</b>	<b>2017-10</b>
Intended use	Anti-Salmonella H:1 records all <i>Salmonella</i> strains, as long as they contain antigen H:1,2, H:1,5, H:1,6, H:1,7. in slide agglutination.
Specificity	No agglutination with heterologous <i>Salmonella</i> H antigens and <i>Salmonella</i> O antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin, 21.11.14  
\_\_\_\_\_  
place, date

  
\_\_\_\_\_  
signature Head of Quality Control