
Certificate of Quality Control

Manufacturer	SIFIN Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany
Product	Enteroclon Anti-Salmonella H z₁₃
Catalogue No. (Package size)	TR 1439 (1 ml)
Lot No.	115 08 12
Expiry date	2015-08
Intended use	Enteroclon Anti-Salmonella H z ₁₃ records all Salmonella strains, as long as they contain antigen H z ₁₃ of the L complex in slide agglutination. The product is not intended for use <i>in vitro</i> for the examination of samples derived from the human body.
Specificity	No agglutination with Salmonella O antigens and heterologous Salmonella H antigens, with the exception of H 1, H z₂₈ and H z₃₉.
Preservatives	Sodium azide 0.9 mg/ ml
Storage	2...8 °C

Berlin,

10.08.12


Head of Quality Control