


## Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN  
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

<b>Product</b>	<b>Anti-Salmonella H:z<sub>13</sub></b>
Catalogue no.	TR1439
Package size	1 ml
<b>Lot no.</b>	<b>8181020</b>
<b>Expiry date</b>	<b>2023-09</b>
Intended use	Anti-Salmonella H:z <sub>13</sub> records all <i>Salmonella</i> strains, as long as they contain antigen H:z <sub>13</sub> of the L complex in slide agglutination. This product is not intended for use <i>in vitro</i> for the examination of samples derived from the human body.
Specificity	No agglutination with <i>Salmonella</i> O antigens and heterologous <i>Salmonella</i> H antigens, <b>with the exception of H:1, H:z<sub>28</sub> and H:z<sub>39</sub>.</b>
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C... 8 °C

Berlin, 04.11.2020  
\_\_\_\_\_  
place, date

  
\_\_\_\_\_  
signature Head of Quality Control