


### Certificate of Quality Control

Manufacturer	<b>SIFIN</b> Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany
Product	<b>Enteroclon Anti-Salmonella H z<sub>38</sub></b>
Catalogue No. (Package size)	TR 1447 (1 ml)
Lot No.	<b>160 12 11</b>
Expiry date	<b>2015-12</b>
Intended use	Enteroclon Anti-Salmonella H z <sub>38</sub> records all Salmonella strains, as long as they contain antigen H z <sub>38</sub> .
Specificity	No agglutination with heterologous Salmonella H antigens and Salmonella O antigens.
Preservatives	Sodium azide 0.9 mg/ ml
Storage	2...8 °C

Berlin, 17.04.12

  
Head of Quality Control