

## Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN  
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485


sifin diagnostics gmbh

Berliner Allee 317-321  
13088 Berlin, Germany

T: +49 30 927 030-0  
F: +49 30 927 030-30  
info@sifin.de  
www.sifin.de

<b>Product</b>	<b>Anti-Shigella sonnei S –form (phase I)</b>
Catalogue no.	TR 2012
Package size	1 ml
<b>Lot no.</b>	<b>97 04 14</b>
<b>Expiry date</b>	<b>2015-04</b>
Intended use	Anti-Shigella sonnei S-form (phase I) records the strains of <i>Shigella sonnei</i> S- form (phase I) in slide agglutination.
Specificity	No agglutination with antigens of Shigella sonnei F-form (phase II), <i>Shigella flexneri</i> , <i>Shigella dysenteriae</i> (type 1 to 10) and <i>Shigella boydii</i> (type 1 to 15).
Preservatives	Sodium azide 0.9 mg / ml
Storage	2 °C...8 °C

Berlin, 08.05.14  
\_\_\_\_\_  
place, date

  
\_\_\_\_\_  
signature Head of Quality Control