


### Certificate of Quality Control

Manufacturer	<b>SIFIN</b> Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany
Product	<b>Enteroclon Anti-Shigella sonnei F-form (Phase II)</b>
Catalogue No. (Package size)	TR 2013 (1 ml)
Lot No.	<b>63 04 12</b>
Expiry date	<b>2015-04</b>
Intended use	Enteroclon Anti-Shigella sonnei F-form (Phase II) agglutinates the F-form (Phase II) of Shigella sonnei
Specificity	No agglutination with Shigella sonnei S-form (Phase I), antigens of Shigella flexneri, Shigella dysenteriae (type 1 to 10) and Shigella boydii (type 1 to 15).
Preservatives	Sodium azide 0.9 mg / ml
Storage	2...8 °C

Berlin,

20.06.12

  
Head of Quality Control