

## Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN  
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh


Berliner Allee 317-321  
13088 Berlin, Germany

T: +49 30 927 030-0  
F: +49 30 927 030-30  
info@sifin.de  
www.sifin.de

<b>Product</b>	<b>Anti-Shigella sonnei S- and F-form (phase I and II)</b>
Catalogue no.	TR 2014
Package size	1 ml
<b>Lot no.</b>	<b>207 10 16</b>
<b>Expiry date</b>	<b>2021-10</b>
Intended use	Anti-Shigella sonnei S- and F-form (phase I and II) records the strains of <i>Shigella sonnei</i> S- (phase I) and F-form (phase II) in slide agglutination.
Specificity	No agglutination with antigens of <i>Shigella flexneri</i> , <i>Shigella dysenteriae</i> (type 1 to 10) and <i>Shigella boydii</i> (type 1 to 15).
Preservatives	Sodium azide 0.9 mg / ml
Storage	2 °C...8 °C

Berlin,  
place, date

4.05.2017

  
signature Head of Quality Control