


Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

Product	Anti-Shigella sonnei S- and F-form (phase I and II)
Catalogue no.	TR2014
Package size	1 ml
Lot no.	5190620
Expiry date	2025-06
Intended use	Anti-Shigella sonnei S- and F-form (phase I and II) records the strains of <i>Shigella sonnei</i> S- (phase I) and F-form (phase II) in slide agglutination.
Specificity	No agglutination with antigens of <i>Shigella flexneri</i> , <i>Shigella dysenteriae</i> (type 1 to 10) and <i>Shigella boydii</i> (type 1 to 15).
Preservatives	Sodium azide 0.9 mg / ml
Storage	2 °C...8 °C

Berlin, 13.07.2020
place, date


signature Head of Quality Control