

### Certificate of Quality Control

Manufacturer	<b>SIFIN</b> Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany
Product	<b>Test Reagent Anti-Coli II</b>
Catalogue No. (Package size)	TR 2121 (1 ml)
Lot No.	<b>27 02 12</b>
Expiry date	<b>2017-02</b>
Intended use	Test Reagent Anti-Coli II records all E. coli strains, as long as they contain the antigens O55:(K59), O86:(K61), O91:(K-), O111:(K58), O119:(K69), O126:(K71), O127:(K63) and O128:(K67).
Specificity	No agglutination with E. coli strains which are recorded by Test Serum /Test Reagent Anti-Coli I and Anti-Coli III.
Preservatives	Sodium azide 0.9 mg/ ml
Storage	2...8 °C

Berlin, 14.12.12

  
\_\_\_\_\_  
Head of Quality Control