

# sifin

## Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN  
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485


sifin diagnostics gmbh

Berliner Allee 317-321  
13088 Berlin, Germany

T: +49 30 927 030-0  
F: +49 30 927 030-30  
info@sifin.de  
www.sifin.de

<b>Product</b>	<b>Anti-Salmonella H:g,m</b>
Catalogue no.	TR 5408
Package size	5 ml
<b>Lot no.</b>	<b>91 04 14</b>
<b>Expiry date</b>	<b>2016-04</b>
Intended use	Anti-Salmonella H:g,m records all <i>Salmonella</i> strains, as long as they contain in an antigen combination the H:g and/or H:m antigens (e. g. H:f,g; H:f,g,s; H:f,g,t; H:g,p; H:g,p,u; H:g,s,t; H:g,t; H:g,z <sub>51</sub> or H:g,m; H:g,m,s,t; H:m,t) in slide agglutination.
Specificity	No agglutination with heterologous <i>Salmonella</i> H antigens and <i>Salmonella</i> O antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin, 10.06.14  
place, date

  
signature Head of Quality Control