

Certificate of Quality Control

Manufacturer	SIFIN Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany
Product	Test serum Anti-Salmonella O 18
Catalogue No. (Package size)	TS 1330 (1 ml)
Lot No.	61 05 11
Expiry date	2016-04
Intended use	Enteroclon Anti-Salmonella O 18 records all Salmonella strains, as long as they contain antigen O 18.
Specificity	No other agglutinations with heterologous Salmonella O antigens and Salmonella H antigens are known.
Preservatives	Sodium azide 0.9 mg / ml
Storage	2...8 °C

Berlin, 08.06.11


Head of Quality Control