

## Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN  
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321  
13088 Berlin, Germany

T: +49 30 927 030-0  
F: +49 30 927 030-30  
info@sifin.de  
www.sifin.de

<b>Product</b>	<b>Anti-Salmonella O:22</b>
Catalogue no.	TS 1332
Package size	1 ml
<b>Lot no.</b>	<b>185 10 14</b>
<b>Expiry date</b>	<b>2016-06</b>
Intended use	Anti-Salmonella O:22 records all Salmonella strains of group G (O:13), as long as they contain antigen O:22.
Specificity	No agglutination with heterologous Salmonella O antigens and Salmonella H antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin, 28.10.14  
place, date

  
signature Head of Quality Control