

sifin

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485


sifin diagnostics gmbh

Berliner Allee 317-321
13088 Berlin, Germany

T: +49 30 927 030-0
F: +49 30 927 030-30
info@sifin.de
www.sifin.de

Product	Anti-Salmonella H:p
Catalogue no.	TS 1414
Package size	1 ml
Lot no.	75 03 14
Expiry date	2019-03
Intended use	Anti-Salmonella H:p records all <i>Salmonella</i> strains, as long as they contain antigen H:p (e. g. H:g,p; H:g,p,s; H:g,p,u; H:m,p,t,[u]) in slide agglutination.
Specificity	No agglutination with heterologous <i>Salmonella</i> H antigens and <i>Salmonella</i> O antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C... 8 °C

Berlin, 10.06.14
place, date


signature Head of Quality Control