

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
 ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321
 13088 Berlin, Germany

T: +49 30 927 030-0
 F: +49 30 927 030-30
 info@sifin.de
 www.sifin.de

Product	Anti-Salmonella H:p
Catalogue no.	TS 1414
Package size	1 ml
Lot no.	102 05 16
Expiry date	2021-05
Intended use	Anti-Salmonella H:p records all <i>Salmonella</i> strains, as long as they contain antigen H:p (e. g. H:g,p; H:g,p,s; H:g,p,u; H:m,p,t,[u]) in slide agglutination.
Specificity	No agglutination with heterologous <i>Salmonella</i> H antigens and <i>Salmonella</i> O antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin,
 place, date

16.08.2016

signature Head of Quality Control

