


Certificate of Quality Control

Manufacturer	SIFIN Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany
Product	Test Serum Anti-Salmonella H q
Catalogue No. (Package size)	TS 1415 (1 ml)
Lot No.	33 03 11
Expiry date	2016-03
Intended use	Enteroclon Anti-Salmonella H q records all Salmonella strains, as long as they contain in an antigen combination the H q antigen (e.g. g,q; g,m,q) in slide agglutination.
Specificity	No agglutination with heterologous Salmonella H antigens and Salmonella O antigens.
Preservatives	Sodium azide 0.9 mg/ ml
Storage	2...8 °C

Berlin, 19.05.11


Head of Quality Control