

### Certificate of Quality Control

Manufacturer **SIFIN** Institut für Immunpräparate und Nährmedien GmbH Berlin,  
Berliner Allee 317-321, 13088 Berlin, Germany

Product **Test Serum Anti-Salmonella H s**

Catalogue No. (Package size) TS 1417 (1 ml)

Lot No. **190 11 13**

Expiry date **2018-07**

Intended use Test serum Anti-Salmonella H s records all Salmonella strains, as long as they contain in an antigen combination the H s (e.g. g,m,s; g,s,t; g,m,s,t; g,s; f,g,s; g,p,s) in slide agglutination.

Specificity no agglutination with Salmonella O antigens and heterologous Salmonella H antigens

Preservatives Sodium azide 0.9 mg/ ml

Storage 2...8 °C

Berlin,

28.11.13

  
Head of Quality Control