

sifin

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh


Berliner Allee 317-321
13088 Berlin, Germany

T: +49 30 927 030-0
F: +49 30 927 030-30
info@sifin.de
www.sifin.de

Product	Anti-Salmonella H:s
Catalogue no.	TS 1417
Package size	1 ml
Lot no.	98 07 13L
Expiry date	2018-07
Intended use	Anti-Salmonella H:s records all <i>Salmonella</i> strains, as long as they contain in an antigen combination the H:s (e. g. g,m,s; g,s,t; g,m,s,t; g,s; f,g,s; g,p,s) in slide agglutination.
Specificity	No agglutination with heterologous <i>Salmonella</i> H antigens and <i>Salmonella</i> O antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin, 06.05.14

place, date



signature Head of Quality Control