

## Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN  
 ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321  
 13088 Berlin, Germany

T: +49 30 927 030-0  
 F: +49 30 927 030-30  
 info@sifin.de  
 www.sifin.de

<b>Product</b>	<b>Anti-Salmonella H:s</b>
Catalogue no.	TS 1417
Package size	1 ml
<b>Lot no.</b>	<b>2570318</b>
<b>Expiry date</b>	<b>2023-03</b>
Intended use	Anti-Salmonella H:s records all <i>Salmonella</i> strains, as long as they contain in an antigen combination the H:s (e. g. g,m,s; g,s,t; g,m,s,t; g,s; f,g,s; g,p,s) in slide agglutination.
Specificity	No agglutination with heterologous <i>Salmonella</i> H antigens and <i>Salmonella</i> O antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin,  
 place, date

13.04.2018

signature Head of Quality Control

