

Certificate of Quality Control

Manufacturer **SIFIN** Institut für Immunpräparate und Nährmedien GmbH Berlin,
Berliner Allee 317-321, 13088 Berlin, Germany

Product **Test serum Anti-Salmonella H t**

Catalogue No. (Package size) TS 1418 (1 ml)

Lot No. **72 08 09**

Expiry date **2014-08**

Intended use Test serum Anti-Salmonella H t records all Salmonella strains, as long as they contain in an antigen combination the H t antigen (e.g. m,t; g,t; g,s,t; f,g,t; g,m,s,t) in slide agglutination.

Specificity No agglutination with Salmonella O antigens and heterologous Salmonella H antigens

Preservatives Sodium azide 0.9 mg/ ml

Storage 2...8 °C

Berlin, 13.10.05


Head of Quality Control