

Certificate of Quality Control

sifin diagnostics gmbh

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

Berliner Allee 317-321
13088 Berlin, Germany

T: +49 30 927 030-0
F: +49 30 927 030-30
info@sifin.de
www.sifin.de

Product	Anti-Salmonella H:t
Catalogue no.	TS 1418
Package size	1 ml
Lot no.	0851217
Expiry date	2022-12
Intended use	Anti-Salmonella H:t records all <i>Salmonella</i> strains, as long as they contain in an antigen combination the H:t (e. g. m,t; g,t; g,s,t; f,g,t; g,m,s,t) in slide agglutination.
Specificity	No agglutination with heterologous <i>Salmonella</i> H antigens and <i>Salmonella</i> O antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin, 25.01.2018
place, date


signature Head of Quality Control