

Certificate of Quality Control

| | |
|---------------------------------|--|
| Manufacturer | SIFIN Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany |
| Product | Test Serum Anti-Salmonella H v |
| Catalogue No. (Package size) | TS 1420 (1 ml) |
| Lot No. | 147 10 12 |
| Expiry date | 2017-10 |
| Intended use | Enteroclon Anti-Salmonella H v records all Salmonella strains, as long as they contain antigen H v (in the antigen combination H I,v). |
| Specificity | No agglutination with Salmonella O antigens and heterologous Salmonella H antigens. |
| Preservatives | Sodium azide 0.9 mg/ ml |
| Storage | 2...8 °C |

Berlin, 23.10.12


Head of Quality Control