

## Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN  
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321  
13088 Berlin, Germany

T: +49 30 700 144-0  
F: +49 30 700 144-30  
info@sifin.de  
www.sifin.de

<b>Product</b>	<b>Anti-Salmonella H:v</b>
Catalogue no.	TS1420
Package size	1 ml
<b>Lot no.</b>	<b>5590720</b>
<b>Expiry date</b>	<b>2025-06</b>
Intended use	Anti-Salmonella H:v records all <i>Salmonella</i> strains, as long as they contain the antigen H:v in the antigen combination H:l,v in slide agglutination.
Specificity	No agglutination with heterologous <i>Salmonella</i> H antigens and <i>Salmonella</i> O antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin,

10.09.2020

place, date

*F. Lenz*

signature Head of Quality Control