

# sifin

## Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN  
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485


sifin diagnostics gmbh

Berliner Allee 317-321  
13088 Berlin, Germany

T: +49 30 927 030-0  
F: +49 30 927 030-30  
info@sifin.de  
www.sifin.de

<b>Product</b>	<b>Anti-Salmonella H:z<sub>6</sub></b>
Catalogue no.	TS 1426
Package size	1 ml
<b>Lot no.</b>	<b>09 01 15</b>
<b>Expiry date</b>	<b>2020-01</b>
Intended use	Anti-Salmonella H:z <sub>6</sub> records all <i>Salmonella</i> strains, as long as they contain antigen H:z <sub>6</sub> in slide agglutination.
Specificity	No agglutination with heterologous <i>Salmonella</i> H antigens and <i>Salmonella</i> O antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin, 26.03.15  
\_\_\_\_\_  
place, date

  
\_\_\_\_\_  
signature Head of Quality Control