

### Certificate of Quality Control

Manufacturer	<b>SIFIN</b> Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany
Product	<b>Test Serum Anti-Salmonella H z<sub>15</sub></b>
Catalogue No. (Package size)	TS 1428 (1 ml)
Lot No.	<b>39 03 13</b>
Expiry date	<b>2018-03</b>
Intended use	Test Serum Anti-Salmonella H z <sub>15</sub> records all Salmonella strains, as long as they contain in an antigen combination the antigen H z <sub>15</sub> (e.g. e,n,x, z <sub>15</sub> ; e,n,z <sub>15</sub> ) in slide agglutination.
Specificity	No agglutination with heterologous Salmonella H antigens and Salmonella O antigens.
Preservatives	Sodium azide 0.9 mg/ ml
Storage	2...8 °C

Berlin,

02.05.13

  
Head of Quality Control