

Certificate of Quality Control

Manufacturer	SIFIN Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany
Product	Test Serum Anti-Salmonella H z₂₄
Catalogue No. (Package size)	TS 1429 (1 ml)
Lot No.	171 11 12
Expiry date	2017-11
Intended use	Test serum Anti-Salmonella H z ₂₄ records all Salmonella strains, as long as they contain antigens H z _{4,Z24}
Specificity	No agglutination with Salmonella O antigens and heterologous Salmonella H antigens.
Preservatives	Sodium azide 0.9 mg/ ml
Storage	2...8 °C

Berlin, 10.12.12



Head of Quality Control