

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
 ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485


sifin diagnostics gmbh

Berliner Allee 317-321
 13088 Berlin, Germany

T: +49 30 927 030-0
 F: +49 30 927 030-30
 info@sifin.de
 www.sifin.de

| | |
|--------------------|---|
| Product | Anti-Salmonella H:z₂₄ |
| Catalogue no. | TS 1429 |
| Package size | 1 ml |
| Lot no. | 134 06 15 |
| Expiry date | 2020-06 |
| Intended use | Anti-Salmonella H:z ₄ ,z ₂₄ records all <i>Salmonella</i> strains, as long as they contain antigens H:z ₄ ,z ₂₄ in slide agglutination. |
| Specificity | No agglutination with heterologous <i>Salmonella</i> H antigens and <i>Salmonella</i> O antigens. |
| Preservatives | Sodium azide 0.9 mg/ml |
| Storage | 2 °C... 8 °C |

Berlin, 29.09.15
 place, date


 signature Head of Quality Control