

Certificate of Quality Control

| | |
|---------------------------------|---|
| Manufacturer | SIFIN Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany |
| Product | Test Serum Anti-Salmonella H 7 |
| Catalogue No. (Package size) | TS 1436 (1 ml) |
| Lot No. | 103 08 11 |
| Expiry date | 2016-08 |
| Intended use | Enteroclon Anti-Salmonella H 7 records all Salmonella strains, as long as they contain antigen H 1,7. |
| Specificity | No agglutination with heterologous Salmonella H antigens and Salmonella O antigens. |
| Preservatives | Sodium azide 0.9 mg/ml |
| Storage | 2...8 °C |

Berlin, 30.09.11



Head of Quality Control