
Certificate of Quality Control

Manufacturer **SIFIN** Institut für Immunpräparate und Nährmedien GmbH Berlin,
Berliner Allee 317-321, 13088 Berlin, Germany

Product **Test Serum Anti-Salmonella H z₂₈**

Catalogue No.
(Package size) TS 1449 (1 ml)

Lot No. **01 01 14**

Expiry date **2017-01**

Intended use Test serum Anti-Salmonella H z₂₈ records all Salmonella strains, as long as they contain antigen H z₂₈ (H ℓ ,z₂₈; H ℓ ,z₁₃,z₂₈).


Specificity No agglutination with Salmonella O antigens and heterologous Salmonella H antigens.

Preservatives Sodium azide 0.9 mg/ ml

Storage 2...8 °C

Berlin,

29.01.14


Head of Quality Control