

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321
13088 Berlin, Germany

T: +49 30 700 144-0
F: +49 30 700 144-30
info@sifin.de
www.sifin.de

Product	Anti-Salmonella H:z₂₈
Catalogue no.	TS1449
Package size	1 ml
Lot no.	6340820
Expiry date	2023-07
Intended use	Anti-Salmonella H:z ₂₈ records all <i>Salmonella</i> strains, as long as they contain antigen H:z ₂₈ in slide agglutination (e. g. H:l,z ₂₈ ; H:l,z ₁₃ ,z ₂₈).
Specificity	No agglutination with heterologous <i>Salmonella</i> H antigens and <i>Salmonella</i> O antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin,

29.09.2020

place, date



signature Head of Quality Control